

# KAI SAFETY KNIVES

STERILE • STAINLESS STEEL • SINGLE USE




## *New System*



Blade protection allows safety handling and avoids accidental injuries  
Precisely engineered design for easy twist and glide action with both hands  
Sterile and ready to use  
Disposable  
5 pieces per dispenser box

## STAB KNIVES

Cataract Surgery / Scleral Tunnel Incision - Pregroove, Side port  
/ Clear Corneal Incision - Side port




<b>SK-A15</b>	Angle 15°	t=0.1mm	<b>SK-A30</b>	Angle 30°	t=0.1mm
<b>SK-A225</b>	Angle 22.5°	t=0.1mm			

Precisely ground in three bevels for ultra sharp incision and piercing.

## CRESCENT KNIVES

Cataract Surgery / Scleral Tunnel Incision - Scleral Tunnel





<b>SK-CPU22</b>	w=2.2mm	Bevel Up	<b>SK-CPU26</b>	w=2.6mm	Bevel Up

Precisely created by computer driven grinding technology for optimal tunnel or flap.

## CLEAR CORNEAL KNIVES

Cataract Surgery / Clear Corneal Incision - Whole Incision




<b>SK-CCR22</b>	w=2.2mm	Double Bevel	<b>SK-CCR28</b>	w=2.8mm	Double Bevel
<b>SK-CCR23</b>	w=2.3mm	Double Bevel	<b>SK-CCR30</b>	w=3.0mm	Double Bevel
<b>SK-CCR24</b>	w=2.4mm	Double Bevel	<b>SK-CCR32</b>	w=3.2mm	Double Bevel

Chemical etching creates ultra sharp blade, producing a double bevel perfectly suited to clear corneal incisions. Unique blade geometry allows high operability for self-sealing corneal tunnel.

## SLIT KNIVES

Cataract Surgery / Scleral Tunnel Incision - Anterior Chamber Entry



<b>SK-SLU22</b>	w=2.2mm	Bevel Up	<b>SK-SLU28</b>	w=2.8mm	Bevel Up
<b>SK-SLU23</b>	w=2.3mm	Bevel Up	<b>SK-SLU30</b>	w=3.0mm	Bevel Up
<b>SK-SLU24</b>	w=2.4mm	Bevel Up	<b>SK-SLU32</b>	w=3.2mm	Bevel Up

Precisely created by computer driven grinding technology to ensure exact incision.

## MVR KNIFE

Cataract Surgery / Scleral Tunnel Incision - Side port  
/ Clear Corneal Incision - Side port



<b>SK-VKB20G</b>	φ 0.9mm