


KAI Surgical Blades & Scalpels



DISPOSABLE SURGICAL BLADES

Regular type

Features


- Available both stainless and carbon steel blades
- Available various blade types
- Sterile and ready to use
- Disposable



Package

- Each blade in a foil package
- 100 blades in a dispenser box
- 100 boxes per shipping carton

Blade type	Stainless steel	Carbon steel
NO.10	310-B	210-B
NO.11	311-B	211-B
NO.12	312-B	212-B
NO.12B	312B-B	212B-B
NO.14	314-B	214-B
NO.15	315-B	215-B
NO.15C	315C-B	215C-B
NO.20	320-B	220-B
NO.21	321-B	221-B
NO.22	322-B	222-B
NO.23	323-B	223-B
NO.24	324-B	224-B
NO.25	325-B	225-B



Actual size

Coated type

Features


- Coated stainless steel blade
- Excellent cutting performance - reducing friction by coating
- Effective for cutting fatty tissue
- Sterile and ready to use
- Disposable

Package

- Each blade in a foil package
- 100 blades in a dispenser box
- 100 boxes per shipping carton

Blade type	Ref no.
NO.10	S310-B
NO.11	S311-B
NO.12	S312-B
NO.14	S314-B
NO.15	S315-B
NO.20	S320-B
NO.21	S321-B
NO.22	S322-B
NO.23	S323-B
NO.24	S324-B
NO.25	S325-B

Cutting force - Nylon thread cutting (dia. 0.1mm)



* Comparison of KAI blades

Surgical knife handles

Autoclavable-Reusable-Plastic

Features

- Polypropylene plastic handle
- With easy blade remover
- Autoclavable
- Non-sterile

Package

- Each handle in a poly bag
- 10 handles in a box
- 100 boxes per shipping carton

Blade type	Ref no.
NO.10s	P-3
NO.20s	P-4

P-3 Handle type : Style NO.3




P-4 Handle type : Style NO.4



Actual size

DISPOSABLE SURGICAL SCALPELS

Regular type



Actual size

Features

- Stainless steel blade
- With transparent cap
- Sterile and ready to use
- Disposable

Package

- Each scalpel in a sterile pouch
- 20 scalpels in a dispenser box
- 50 boxes per shipping carton



Blade type	Ref no.
NO.10	510-A
NO.11	511-A
NO.12	512-A
NO.14	514-A
NO.15	515-A
NO.20	520-A
NO.21	521-A
NO.22	522-A
NO.23	523-A
NO.24	524-A
NO.25	525-A

Coated type

Features


- Coated stainless steel blade
- Excellent cutting performance - reducing friction by coating
- Effective for cutting fatty tissue
- Sterile and ready to use
- Disposable

Package

- Each scalpel in a sterile pouch
- 20 scalpels in a dispenser box
- 50 boxes per shipping carton

Blade type	Ref no.
NO.10	S510-A
NO.11	S511-A
NO.12	S512-A
NO.14	S514-A
NO.15	S515-A
NO.20	S520-A
NO.21	S521-A
NO.22	S522-A
NO.23	S523-A
NO.24	S524-A
NO.25	S525-A

Safety type



Actual size

Features

- Safe from accidental needle sticks
- Blade extends and retracts with one hand
- Glides easily and comfortably into locked position
- Sterile and ready to use
- Disposable

Package

- Each scalpel in a sterile pouch
- 10 scalpels in a dispenser box
- 100 boxes per shipping carton



Blade type	Ref no.
NO.10	SS-10
NO.11	SS-11
NO.12	SS-12
NO.14	SS-14
NO.15	SS-15
NO.20	SS-20
NO.21	SS-21
NO.22	SS-22
NO.23	SS-23
NO.24	SS-24
NO.25	SS-25



kai industries co., ltd.

1110, Oyana, Seki City, Gifu Pref., 501-3992, Japan
Phone: +81 (0) 575 28 5040
Fax: +81 (0) 575 28 6611
E-mail: info-ind@kai-group.com
www.kai-group.com/global/en/



kai Europe GmbH

Kottendorfer Straße 5, 42697 Solingen, Germany
Phone: +49 (0) 212 23 23 80
Fax: +49 (0) 212 23 23 899
E-mail: info@kai-europe.com
www.kai-europe.com

Product specifications may change without notice. ■